



SYDMED
A U S T R A L I A



NuWav[™]


360° Endo-Cavity 12MHz Ultrasound Probe for Anal / Sphincter Trauma

Laborie offers the NuWav™ 360° Endo-Cavity 12MHz Ultrasound Probe, which allows imaging of key structures to diagnose anal sphincter competence.



Features:

- Easy to use - Simply plugs into your laptop's USB port. Images are displayed on the laptop's screen
- Smallest Transrectal Ultrasound (TRUS) System available worldwide!
- Affordable imaging diagnostics for anal sphincter competence
- Screening procedures looking at muscle bulk, symmetry, and potential tears/damage made easy and affordable
- Perfect adjunct to Anorectal Manometry (ARM) diagnostics available through Laborie's global leading UDS® software platform
- Self-contained hardware allows for simple cleaning, service, and upgrade
- Software controls for TGC, Brightness and Contrast, Depth, Power and Frequency
- Variable frequencies of 5MHz, 7.5MHz, 12MHz



Ultrasound image showing tear of IAS and EAS

SYDMED

Suite 301/30-40 Harcourt Parade

Rosebery NSW 2018 Australia

Tel: (02)9317 4777 Fax: (02) 9317 3165

Email: info@sydmed.com

Web: <http://www.sydmed.com>