



Gynecology Procedure Chair

SONESTA **6302**



The 6302 plus...
dual-motor, equals
affordability and
versatility for just
about any type
of Gynecology
procedure.



STILLE



Gynecology Procedure Chair

SONESTA 6302

There Is A 6302 That Fits Your Needs.

- Urodynamics
- General
- Examinations
- Cystoscopy
- Ultrasound
- Biofeedback
- Vasectomy
- OB/GYN Procedures



Engineered To Be Flexible.

- Stille flexibility for office, clinic, or hospital
- Easy positioning for patient and physician comfort
- Accepts the full range of Stille Accessories
- Wide range of adjustment, from a low of 16" (40cm) to 39.4" (100cm) high
- Lightweight and mobile easy to move to other rooms; unmatched access for floor cleaning and disinfecting.
- Tough and reliable supports up to 375 lbs - 170 kg for years of service in the most demanding practices.

The 6302 is engineered to be the affordable, flexible solution in your examination room.

This two-motor table offers complete flexibility for offices on a budget. The motorized height & seat positioning lets you comfortably position your patients while the gas piston assist controls the associated back angle. The 6302 includes a list of standard accessories that covers many day-to-day procedures or can be customized for more specific needs from a full line of Stille upgrade accessories.

SPECIFICATIONS:

LENGTH WITH EXTENSION CUSHION	175cm	69"
LENGTH WITHOUT	125cm	49.2"
WIDTH	65cm	25.5"
HEIGHT (MIN)	40cm	16"
HEIGHT (MAX)	100cm	39"
LIFT CAPACITY	170kg	375lbs
WHEEL DIAMETER (LOCKING)	10cm	4"
NET WEIGHT	<100kg	<220lbs

ELECTRICAL:

120/230v 24v DC Output
Meets EMC (Electromagnetic Compatibility) for IEC 60601
All Products carry the CE mark.

Standard Accessories:

2 motor foot control	524- 6353
Cushion w/170mm cutout	524-6366-3
Paper roll holder	524-6326
2 heavy duty Rissler leg rests	524-6336
2 leg rest holders	524-6338
Drainage Receptacle	524-6325-1
2 hand grips/arm rests	524-6339-2

DISTRIBUTED IN U.S. BY:

Laborie Medical

400 Avenue D, Suite 10
Williston, VT 05495
Phone: +1 800.522.6743
Fax: +1 802.878.1122
Web: www.laborie.com

CONTACT IN CANADA:

Laborie Medical

2101 Boulevard Lapinière
Brossard, Québec J4W 1L7
Phone: +1 888.522.6743
Fax: +1 450.671.7182
Web: www.laborie.com

